

RICHMED®

TAIZHOU RICH MEDICAL PRODUCT CO.,LTD



RICH MEDICAL PRODUCTS MANUFACTURER



0197

ISO 13485:2016

ISO 11135 2014

◆ Customer's satisfaction is always the company's first priority.

Taizhou Rich Medical Products Co., Ltd, founded in 2003, specialize in the manufacture and trade of disposable medical products. Our company covers an area of 13,000 square meters. Among it there are nearly 2,000 square meters of 100,000 level dust-free workshops and an ethylene oxide sterilization workshop. Our production lines consist of four major categories: medical face masks, medical dressing kits, absorbent cotton products, and medical plastic products. Company currently has CE and FDA certificates. Also we have QMS certificates of ISO13485 and ISO11135. At the same time, we also have the medical device production license and medical products registration certificates in China. The pursuit of "integrity, quality, and innovation" is our common goal.

Let RICHMED and employees grow together, we will continue to keep close cooperation with customers in wider areas, to supply more high-quality medical products, benefit the human health.





QUALITY SAFETY COMFORT



RICH MEDICAL

PART ONE ABSORBENT COTTON

Cotton balls.....	1
Absorbent cotton wools	3
Cotton tipped applicators.....	5
Cotton zig zag.....	7
Cotton sliver(cotton coil)	8
Gauze.....	9
Non-wove sponges	10
Sterile non-wove sponges	11
Cosmetic cotton pad	12
Eye pad	13
Gamgee roll/pad	14
Medical Underpads(Soft Care pad)	15



PART TWO MEDICAL DRESSING KITS

Basic Dressing Set	17
Ophthalmic Dressing Pack.....	18
Circumcision Dressing Pack.....	19
Disposable oral care dressing pack.....	20
IV Start Kit	20
Vaginal Examination Set	22
Suture Removal Kit.....	23





CONTENTS

PART THREE MEDICAL PLASTIC

Sponge cleaning swabs	24
Disposable plastic forceps	24
Disposable plastic sponge	25
Other plastic consumables	26
Medical Plastic tray	28
Surgical hand brush	29
CPR Breathing Mask	30



PART FOUR MEDICAL/SURGICAL FACE MASK

Medical/Surgical Face Mask	31
Mask Production	32
Mask detection machine	33
Ethylene Oxide Sterilizer	34



PART FIVE FIRST AID KIT

First Aid Kit	35
---------------	----



Cotton Balls



ABSORBENT
COTTON BALL
Color cotton ball





name	spec	packing	dimension (mm)
Cotton ball	0.3g/ball	50 g/bag 12bags/ box 8boxes/ ctn	560x400x550
	0.5g/ball	500 g/bag 20bags/ ctn	640x460x480
	1.0g/ball	100g/bag 50bag/ ctn	550x380x430
	2.0g/ball	500balls/bag 5bags/ctn	530x380x420
	3.0g/ball	400balls/bag 5bags/ctn	530x380x500
	4.0g/ball	250balls/bag 5bags/ctn	530x380x420
	5.0g/ball	200balls/bag 5bags/ctn	530x380x420
	10g/ball	200balls/bag 4bags/ctn	640x460x400
Color cotton ball	0.5g/ball	50g/bag 200bags/ ctn	640x460x480
	1.0g/ball	100g/bag 100bags/ ctn	640x460x480
Sterile cotton ball	0.3g/ball	10 balls/pack 100packs/ bag 10bags/ ctn	580x400x450
	0.5g/ball	5 balls/pack 1000packs/ ctn	530x400x450
	0.7g/ball	10 balls/pack 1000bags/ ctn	730x460x460
	1.0g/ball	5 balls/pack 20packs/bag 50bags/ ctn	540x460x480
	5.0g/ball	15 balls/bag 25bags/ctn	430x320x380



◆ Absorbent Cotton Wools



◆ Dental Cotton Roll



spec	packing			dimension(mm)
8x38mm	50pcs/bundle	20bundle/bag	20bags/ctn	370 x 230 x370
10x38 mm	50pcs/bundle	20bundle/bag	20bags/ctn	560 x 420 x275
12x38mm	50pcs/bundle	20bundle/bag	20bags/ctn	540 x 350 x460



◆ Cotton Rolls



spec	packing	dimension(mm)
25 gram	200 roll/ctn	370 x 230 x370
50 gram	200 roll/ctn	560 x 420 x275
100 gram	100 roll/ctn	540 x 350 x460
200 gram	40 roll/ctn	460 x 400 x320
250 gram	40 roll/ctn	460 x 350 x360
400 gram	25 roll/ctn	420 x 310 x420
454 gram	25 roll/ctn	450 x 310 x440
500 gram	25 roll/ctn	470 x 310 x460
1000 gram	12 roll/ctn	600 x 310 x440
4000 gram	6 roll/ctn	910 x 310 x450

◆ Cotton Tipped Applicators



◆ Wooden Stick





◆ Plastic Stick



name	spec (single head)	packing			dimension (mm)
Wooden stick	3" φ 5mm	100pcs/bag	10bags/box	10boxes/ctn	420x200×250
	6" φ 5mm	100pcs/bag	10bags/box	20boxes/ctn	430x370×250
	6" φ 10mm	50pcs/bag	20bags/box	8 boxes/ctn	630x370×490
	8" φ 10mm	50pcs/bag	20bags/box	8 boxes/ctn	570x400×480
	8" φ 20mm	20pcs/bag	20bags/box	8 boxes/ctn	540x400×480
Plastic stick	3" φ 5mm	100pcs/bag	10bags/box	10 boxes/ctn	420x200×250
	6" φ 5mm	100pcs/bag	10bags/box	20 boxes/ctn	430x370×370
	6" φ 10mm	50pcs/bag	20bags/box	8 boxes/ctn	630x370×490
	8" φ 10mm	50pcs/bag	100bags/ctn		520x500×320
	8" φ 20mm	20pcs/bag	20bags/box	8 boxes/ctn	570x400×480
Sterile cotton applicator	3" φ 5mm	5pcs/pack	100bags/bag	10 bags/ctn	400x135×300
	6" φ 5mm	5pcs/pack	100bags/bag	10 bags/ctn	420x370×230
	6" φ 10mm	2pcs/pack	25packs/bag	45 bags/ctn	580x430×415
	8" φ 10mm	1pc/bag	100bags/box	50 boxes/ctn	630x500×330
	8" φ 20mm	1pc/bag	100bags/box	50 boxes/ctn	730x500×630

◆ Cotton Zig Zag



spec	packing	dimension(mm)
5 gram	100 bags/ctn	260 x 130 x160
10 gram	100 bags/ctn	310 x 260 x230
20 gram	100 bags/ctn	460 x 310 x260
25 gram	100 bags/ctn	490 x 200 x330
50 gram	100 bags/ctn	470 x 320 x460
100 gram	100 bags/ctn	490 x 470 x560
200 gram	60 bags/ctn	480 x 300 x490
250 gram	40 bags/ctn	500 x 480 x460
500 gram	10 bags/ctn	510 x 360 x320
1000 gram	5 bags/ctn	580 x 510 x180



◆ Cotton Sliver(Cotton Coil)



name	spec	packing	dimension (mm)
Cotton Coil	2-9g/m	500g/bag 10bags/ctn	540x520×350
		10kg/bags	Dia:500 H:930

Gauze



name	spec	packing	dimension (mm)
Gauze Swab	5×5cm 8ply 13thread	100pcs/bag 200bags/ctn	530x530×140
	7.5×7.5cm 8ply 13thread	100pcs/bag 100bags/ctn	390x390×260
	10×10cm 8ply 13thread	100pcs/bag 50bags/ctn	570x520×140
Gauze roll	36 " ×100y 13thread	20rolls/ctn	750x320×360
	36 " ×100m 13thread	20rolls/ctn	790x340×360
Gauze bandage	5cm×5m 13thread	100pcs/ctn	530x260×220
	7.5cm×5m 13thread	100pcs/ctn	770x260×220
	10cm×5m 13thread	100pcs/ctn	580x510×220



name	spec	packing	dimension (mm)
Sterile Gauze Swab	10×10cm 8ply 13thread	5pcs/bag 20pouches/box 20boxes/ctn	570x290×330
	7.5×7.5cm 8ply 13thread	5pcs/bag 20pouches/box 40boxes/ctn	520x420×260



◆ Non-Wove Sponges



name	spec	packing	dimension(mm)
Non Woves Sponges	2"×2" 4ply 30gsm	100pcs/pack 100packs/ctn	280×280×420
	3"×3" 4ply 30gsm	100pcs/pack 100packs/ctn	390×390×420
	3"×3" 6ply 30gsm	100pcs/pack 50packs/ctn	390×390×315
	3"×3" 8ply 30gsm	100pcs/pack 50packs/ctn	390×390×420
	3"×3" 12ply 30gsm	100pcs/pack 50packs/ctn	625×390×390
	4"×4" 4ply 30gsm	100pcs/pack 100packs/ctn	550×550×420
	4"×4" 6ply 30gsm	100pcs/pack 50packs/ctn	550×550×320
	4"×4" 8ply 30gsm	100pcs/pack 50packs/ctn	550×550×420
	4"×8" 10pcs/pack	40packs/ctn	300×450×230
	6"×9" 5pcs/pack	80packs/ctn	430×230×280
	6"×18" 5pcs/pack	80packs/ctn	430×250×320
	12"×12" 5pcs/pack	80packs/ctn	400×210×420
	12"×18" 5pcs/pack	40packs/ctn	350×280×280
	15"×15" 5pcs/pack	40packs/ctn	420×260×280
	18"×18" 5pcs/pack	40packs/ctn	400×400×230



◆ Sterile Non-Woven Sponges



name	spec	packing			dimension(mm)
Sterile Non-Wove Sponges	2"×2" 4ply	5pcs/pack	50packs/box	30boxes/ctn	470×320×500
	3"×3" 4ply	5pcs/pack	50packs/box	30boxes/ctn	570×400×500
	4"×4" 4ply	5pcs/pack	50packs/box	24boxes/ctn	580×530×500



◆ Cosmeric Cotton Pad



name	spec	packing		dimension (mm)
Round pad	φ 58mm	50pcs/bag	48bags/ctn	350x250×240
	φ 58mm	80pcs/bag	48bags/ctn	430x250×350
	φ 58mm	100pcs/bag	48bags/ctn	470x250×350
Square pad	5×5cm	450g/bag 1000/bag	24bags/ctn 10bags/ctn	630x500×480 770x420×410
	8×10cm			
	7.5×7.5cm			
	10×10cm			
	15×15cm			

◆ Eye Pad



name	spec	packing	dimension (mm)
Non-woven	6×8cm	1pc/bag 50bags/box 50boxes/ctn	500x425×360
	5×7cm	1pc/bag 50bags/box 50boxes/ctn	625x500×360
	5.4×6.7m	1pc/bag 50bags/box 50boxes/ctn	625x500×360
	4.1×6.7m	1pc/bag 50bags/box 50boxes/ctn	625x500×360
Gauze	6×8cm	1pc/bag 50bags/box 50boxes/ctn	500x425×360
	5×7cm	1pc/bag 50bags/box 50boxes/ctn	625x500×360
	5.4×6.7cm	1pc/bag 50bags/box 50boxes/ctn	625x500×360
	4.1×6.7cm	1pc/bag 50bags/box 50boxes/ctn	610x420×310



◆ Combine Dressing Roll




Absorbent Gauze Cover



Gamgee



Non-Woven Cover



name	spec	packing	dimension (mm)
Gauze	30cm×3m	12roll/ctn	400x330×300
	45cm×1.5m	12roll/ctn	480x270×340
	45cm×2.7m	12roll/ctn	480x270×340
	40cm×2m	16roll/ctn	420x330×320
	40cm×5m	8roll/ctn	600x415×300
	20cm×5m	16roll/ctn	610x340×420
Non-woven	9cm×10m	30 roll/ctn	490x480×330
	20cm×10m	16 roll/ctn	610x420×310
	30cm×4m	20 roll/ctn	650x320×520
Gamgee	25cm×0.4m	1pc/pack 80packs/ctn	550x450×440



◆ Medical Underpads(Soft Care pad)



name	spec	packing	dimension (mm)
Medical Underpads (Soft Care pad)	60×90cm	10pcs/bag 8bags/ctn	500×430×590



◆ Medical Dressing Kits



◆ MISSION

Make every effort to continuously contribute to the improvement of quality, comfort, performance and safety of care in the hospital, as the consulting room or at home.

◆ Basic Dressing Set



name	introduce
Basic Dressing Set	1 x Dressing tray
	5 x Cotton wool balls middle
	5 x Gauze swab
	3 x Dressing forceps
	1 x Surgical drape
	1 x Hand towel
	1 x Limpet bag
Dimension	530x330x370mm(100 sets)



◆ Ophthalmic Dressing Pack



name	introduce
Ophthalmic Dressing Pack	1 x Blister 18.5 x 15.5 x 4.5
	4 x Non-woven swab 7.5cm x 7.5cm
	1 x Eye pad
	1 x Hand tower
	1 x Sterile field 45 x 45cm
	5 x Cotton ball 0.5 gram/pc
	1 x Red limpt bag with tie
Dimension	530 x 330 x 370 mm(100 kits)

◆ Circumcision Dressing Pack



name	introduce
Circumcision Dressing Pack	1 x Dressing tray 2 compartments 15x18x4.5cm
	8 x Cotton wool balls 0.7g
	5 x Gauze swabs 10cm-8ply 13thread
	1 x Forceps 13cm
	1 x Surgical drape 40cmx40cm
	1 x Pair latex gloves
	1 x Limpet bag red 30x43cm
Dimension	520x415x450mm (100 kits)

MEDICAL DRESSING KITS RICHMED®



Disposable Oral Care Dressing Pack



name		introduce
Disposable oral care dressing pack		1 x Hand towel
		1 x Plastic hemostat
		1 x Plastic fprceps
		1 x Dressing towel
		1 x Tongue depressor
		14 x Cotton balls
		1 x Pouch Lube swab
Dimension		400 x 370 x 360 mm (100 kits)

◆ Iv Start kit



name	introduce
Iv Start Kit	1 x pc Latex Free Tourniquet 3/4" x 18"
	1 x pc Surgical tape 2cm x 1 m
	1 x pc IV Wound dressing 6 x 7 cm
	1 x pc Alcohol swab pad
	1 x pc Povidone -iodine prep pad
	2 x pcs Gauze swab 5cm x 5cm - 8ply
Dimension	430 x 240 x 330 mm (200 kits)



Vaginal Examination Set



name		introduce
Vaginal Examination Set		1 x Plastic Gallipot 11.5 x 7.5 x 4 cm
		1 x pc Hand tower
		1 x pc Plasric Patient Drape
		1 x pc Yellow limpet bag
		1 x pc Large cotton balls
		1 x pc Sanitary pad
		1 x pc Outer Wrap-Sterile Field
Dimension		350 x 340 x 320 mm (30 sets)

◆ Suture Removal Kit



name	introduce
Suture Removal Kit	1 x suture removal scissors
	1 x plastic forceps 10.5cm
	1 x gauze sponge 7.5cm-8ply 13thread
Dimension	330 x 320 x 290 mm (100 kits)



◆ Spong Cleaning Swabs




name	spec	packing		dimension (mm)
Oral cleaning swab	13cm	20pcs/bag	50bags/ctn	380x280x300
	15cm radiused	20pcs/bag	50bags/ctn	560x500x380
	15cm cornered	20pcs/bag	50bags/ctn	570x550x360
Sponge clean stick	18cm radiused	20pcs/bag	50bags/ctn	560x500x380
	18cm cornered	20pcs/bag	50bags/ctn	610x550x430

◆ Disposable Plastic Forceps

name	spec	packing		dimension (mm)
1#	10.8cm	500pcs/bag	4bags/ctn	460x300x320
3#	11.2cm	500pcs/bag	4bags/ctn	460x380x310
4#	11.2cm	500pcs/bag	4bags/ctn	460x380x310
5#	12.5cm	500pcs/bag	4bags/ctn	460x380x310
6#	18cm	1000pcs/ctn		460x380x310
8#	12cm	500pcs/bag	4bags/ctn	460x380x310
9#	13cm	250pcs/bag	4bags/ctn	380x300x320
10#	14cm	250pcs/bag	4bags/ctn	400x320x320
15#	12.5cm	250pcs/bag	4bags/ctn	380x300x320
18#	12.6cm	500pcs/bag	4bags/ctn	460x380x310



◆ Disposable Plastic Sponge



name	spec	packing	dimension (mm)
Haemostatic forceps	14cm	500pcs/ctn	380x280×230
Towel forceps	18cm	500pcs/ctn	380x320×280
	24cm	500pcs/ctn	380x370×280



name	spec	packing	dimension (mm)
Kindey Basin	500ml 22g/pc	250pcs/ctn	230x230×440
	700ml 33g/pc	100pcs/ctn	280x260×250
	700ml 18g/pc	300pcs/ctn	600x400×260



◆ Other Plastic Consumables




name	spec	packing
Medicine cup	60ml	500pcs/ctn
	120ml	
	250ml	
	450ml	
	1000ml	



name	spec	packing
Intestinal irrigator	150ml	individual paper box



name	spec	packing	dimension (mm)
Tick Remover/tweezer	10.9cm	500pcs/bag 4bags/ctn	380x300×310
		1pc/card 12cards/box 12boxes/ctn	340x200×280



name	spec	packing	dimension (mm)
Vomit bag	1000ml	200pcs/ctn	290x210×260



◆ Medical Plastic Tray



Transparent 24#



Whitening 32#



Whitening 33#



Whitening 34#



Transparent 27#



Whitening 37#



Whitening 35#



Whitening 38#

name	spec	packing	dimension (mm)
Transparent 24#	27x19x4cm	200pcs/ctn	400x290x440
Transparent 27#	11.5x8.4x0.6cm	1000pcs/ctn	470x250x450
Whitening 32#	18.2x13.5x14.2cm	500pcs/ctn	560x230x470
Whitening 33#	18.2x13.5x6.3cm	500pcs/ctn	430x365x525
Whitening 34#	22.7x17.7x5cm	270pcs/ctn	560x260x245
Whitening 35#	27.5x22.2x6.1cm	136pcs/ctn	570x290x240
Whitening 37#	11.0x7.8x3.3cm	1000pcs/ctn	460x220x420
Whitening 38#	11.5x11.5x2.3cm	600pcs/ctn	370x230x480

◆ Surgical Hand Brush

RH160004



RH160001



RH160003




RH160002



name	spec	packing dimension (mm)	
Surgical Hand Brush	Sponge Brush 8x5x2.1cm Iodine 20ml	250pcs/ctn	480x380x330
Surgical Hand Brush With Iodine 7.5%/10%	Sponge Brush 8x5x2.1cm Iodine 20ml	300pcs/ctn	580x375x400
Surgical Hand Brush With Chiorhexidine4%	Sponge Chiorhexidine 20ml	300pcs/ctn	580x375x400





◆ CPR Breathing Mask



name	spec	packing	dimension (mm)
CPR Breathing Mask Meltblown	20x30cm	1pc/bag 1000bags/ctn	380×380×310
CPR Breathing respirator	1#20x20cm	1pc/bag 500bags/ctn	370×300×280
	2#20x20cm	1pc/bag 500bags/ctn	510x350×320
	3#20x20cm	1pc/bag 500bags/ctn	360×200×240



Medical/Surgical Face Mask



name	spec	bfe	packing	dimension (mm)
Medical Face Mask	Tie-loop stype Ear-loop stype	≥95%	10pcs/pkt 25pkts/bag 8bags/ctn 50pcs/box 40boxes/ctn	170x100 145x100
Surgical Face Mask	Tie-loop stype Ear-loop stype	≥98%/99%	10pcs/pat 25pkts/bag 8bags/ctn 50pcs/box 40boxes/ctn	175x95 142x90



◆ Mask Production



◆ Mask Detection Machine



PFE



Gas exchange pressure
difference tester



Flame-Resistance Tester




BFE



Synthetic blood penetrometer



Microbiology





TAIZHOU RICH MEDICAL PRODUCTS CO.,LTD

LISA XIA · Marketing Director

Mobile / WhatsApp / WeChat: +86 150 6108 8399

Email: lisa.medrainbow@gmail.com

Website: www.richmedcare.com

Add: No. 20, Zijing Road, Suchen Town, Taizhou, Jiangsu, China